CE  
0197

**KARDIO-TEST**®  
M E D I C A L

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**KARDIO-TEST**<sup>®</sup>  
M E D I C A L

# Instruction Manual

Digital blood pressure  
and heart rate monitor KTB-02

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Established in 1949

**TECH - MED**

®

- Thank you for choosing the KARDIO-TEST MEDICAL KTB-02 digital blood pressure and heart rate monitor.
- Please read these operating instructions carefully. This will ensure safe use of the product. Please keep the operating instructions.

## Important safety precautions

### Read this before use

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or damaged.
4. Do **NOT** use the equipment where aerosol sprays are being used, or where oxygen is being administered.
5. Do **NOT** use under any circumstances on newborns.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
7. Keep the equipment and its flexible cord away from hot surfaces.
8. Do **NOT** apply the cuff to areas other than the place directed.
9. Proper maintenance and periodically calibration are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact local customer service for help.

## BEFORE YOU BEGIN

### 1. Intended use.

The Blood Pressure Monitor is a system designed to measure people's systolic and diastolic blood pressure and pulse rate at home and in clinical settings by using a non-invasive technique in which an inflatable cuff is wrapped on the upper arm.

This system is intended for individuals age 18 or above.

Do not use it for diagnosis of hypertension or testing on newborns, babies, young children or persons who cannot express their consent.

### 2. Test principle.

Blood pressure is measured non-invasively at the arm based on oscillometric method.

If you have been diagnosed with a severe arrhythmia or irregular heartbeat, atrial or ventricular premature beats or atrial fibrillation, measurements made with this instrument should only be evaluated after consultation with the doctor or your healthcare professionals.

#### **Electromagnetic interference**

The device contains sensitive electronic components. Avoid strong electromagnetic fields in the direct vicinity of the device (e.g. mobile phones, microwave ovens). Electrical interference can lead to temporary impairment of the measuring accuracy).

### 3. Monitor overview

#### a) Device


#### b) Cuff

We offer a standard cuff, measuring size: 22-32cm, suitable for most of arm circumferences

Do not use accessories or equipment which are not certified by manufacturer, they might increase emissions or decrease immunity of this system

#### 4. Display screen

- |                 |                                      |
|-----------------|--------------------------------------|
| 1. Battery icon | 9. Time                              |
| 2. Speaker      | 10. Systolic pressure                |
| 3. User icon    | 11. Diastolic pressure               |
| 4. Cuff check   | 12. Unit icon                        |
| 5. Moving check | 13. Heartbeat icon                   |
| 6. WHO bar      | 14. Irregular heartbeat detector     |
| 7. Memory icon  | 15. Heartbeat icon                   |
| 8. Memory index | 16. Average memory of latest 3 times |

## SETTING THE MONITOR

Before using your monitor for the first time, check and update these settings. Make sure you complete the steps below and save the desired settings.

#### 1. Insert the batteries

After you have unpacked your device, first insert the attached batteries (4 x size AAA 1.5V) with the indicated polarity. The battery compartment is located at the bottom of the device (see illustration).

#### Note


- After the warning  appears, device is blocked till the batteries have been replaced.
- Please use 1.5V AAA size Batteries
- If the blood-pressure monitor is left unused for long periods, please remove the batteries.

## 2. Setting mode

### 1. User setting


When LCD is off, press and hold the SET button for 3 seconds,  
The LCD will display user icon, Press the MEM button to select user.


### 2. Setting the date and time, voice and switch unit

After confirm the user, Press the SET button. The LCD will display flashing years. Press MEM button until the correct year appears. Press SET button to confirm and change to Month setting.

b. When the month flashes, press MEM button until the correct month appears Press SET button to confirm and change to date setting.

c. Set the day, hour and minute in the same way as described above

d. After minute setting, The device enter into voice switch setting state,  
Press The MEM button to change the voice status.

e. Press the SET button to confirm the voice switch. the device enter into UNIT selection. Press MEM button to change unit, PA ON is kPa unit, PA OFF is mmHg unit.

### 3. Clear Memories

In the shutdown state, Press and hold the MEM button about 6 seconds, Then press the SET button, The LCD will display ALL EE, All users' memories will be deleted.

Before you delete all the readings stored in the memory, please make sure you will not need to use them any more at a later time

### Note

If the meter is idle for 1 minutes during the setting mode, it will switch off automatically. Press START/ STOP button to make measurement at any time when setting.

## TESTING YOUR BLOOD PRESSURE

### 1. Before measurement

- Avoid caffeine, tea, alcohol and tobacco for at least 30 minutes before measurement.
- Wait 30 minutes after exercising or bathing before measurement Sit or lie down for at least 10 minutes before measurement.
- Do not measure when feeling anxious or tense.
- Take a 5-10 minutes break between measurements. This break can be longer if necessary, depending on your physical condition.
- Keep the records for your doctor as reference.
- Blood pressure naturally varies between each arm. Always measure your blood pressure on the same arm.

### 1. Fitting the cuff properly

- a. Plug the cuff connector to the cuff socket on the device left side
- b. Assemble the cuff as shown below. The smooth surface should be inside of the cuff loop and the metal D-ring should not touch your skin
- c. Stretch your left(right)arm in front of you with your palm facing up. Slide and place the cuff onto your arm to let the air tube and artery mark region toward the lower arm Wrap and tighten the cuff above your elbow. The artery mark on the edge of the cuff should be approx. 0.8" to 1.2" (2 cm to 3 cm) above your elbow. Align the tube over the main arteries on the inside.

- d. Leave a little free space between the arm and the cuff: you should be able to fit two fingers between them  
Remove all clothing covering or constricting the measurement arm
  - e. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm
2. Proper measurement position
- a. Sit down for at least 10 minutes before measuring
  - b. Place your elbow on a flat surface. Relax your hand with the palm facing up
  - c. Make sure the cuff is about the same height as the location of your heart.
  - d. Press the start/stop button to switch on the monitor
  - e. Remain still and do not talk or move during the measurement


### Warning

If the cuff is relatively lower (higher) than the heart, the obtained blood pressure value could be higher (lower) than the actual value. A 15 cm difference in height may result in an error around 10 mmHg.

2. Taking measurements
  1. Taking a measurement

**Always apply the pressure cuff before turning on the monitor.**

- a. Press START/STOP button, pump begins to inflate the cuff. The rising pressure in the cuff is shown on the display.

- b. Pressure in cuff increases to working pressure. Then pump stops and pressure falls. Note: In case the inflation pressure is not sufficient to determine measurement result, the monitor automatically re-inflates to a higher level.
- c. The heart symbol „♥” starts to flash every heart beat once pulse is detected.

- d. After the measurement, systolic pressure, diastolic pressure, pulse rate and color bars classified from ESH (European Society Hypertension) standards shown on the display.
- e. If irregular heart beat was detected during the measurement, LCD displays „⚡” icon to remind users of heartbeat irregularity. Please consult your doctor if it appears often.


**Note:**

The three-color backlight indication function ONLY for some special models

1. Red backlight lasts when the reading of systolic is  $\geq 160$ mmHg or diastolic reading  $\geq 100$  mmHg
2. Yellow backlight lasts when the reading of systolic is higher than 130mmHg but less than 160mmHg, diastolic reading higher than 85mmHg but less than 100mmHg in the meantime
3. Green backlight lasts when the reading of systolic is less than 130mmHg, as well as diastolic reading is less than 85mmHg
4. When check the memory, the backlight colour is the same as the colour of the end of the measurement

### 1. Average of last three memories

Press MEM button, LCD shows a set of values averaged of latest three times memories, accompanied by an icon "AVG"


### 2. Memory recall

The device stores 99 most recent blood pressure test results along with respective dates and times in the monitor memory (each user 99 memories).

Reviewing the results

- Press MEM button again after you see the average of all stored last three memories, the single memory reading appears.
- Per the setting, the first reading you see is the last blood pressure result. Press MEM button again, the next measurement eg. M90, as well as the previous 89 measurements (M90, M89... M1) can be displayed one after the other.

### Note:

The three-color backlight indication function ONLY for some special models

- Red backlight lasts when the reading of systolic is  $\geq 160$  mmHg or diastolic reading  $\geq 100$  mmHg
- Yellow backlight lasts when the reading of systolic is higher than 130 mmHg but less than 160 mmHg, diastolic reading higher than 85 mmHg but less than 100 mmHg in the meantime
- Green backlight lasts when the reading of systolic is less than 130 mmHg, as well as diastolic reading is less than 85 mmHg
- When check the memory, the backlight colour is the same as the colour of the end of the measurement

### 1. Caring for your monitor

To avoid the monitor attracting dirt, dust or other contaminants, wash and dry hands thoroughly before use.

#### Cleaning

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the monitor.
- Do NOT wash the pressure cuff.
- Do NOT iron the pressure cuff.

#### Monitor Storage

- Storage condition: -20°C to 55°C (-4°F to 131°F), below 95% relative humidity.
- Always store or transport the monitor in its original storage case.
- Avoid droppiing or heavy impact.
- Avoid direct sunlight and high humidity.

### 1. Reference values

Clinical studies show that adult diabetes is often accompanied by elevated blood pressure.

People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment

Monitoring your routine blood pressure trend helps you to know your body condition, Human blood pressure naturally increases after reaching middle age This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol(LDL) adhering to the blood vessels, diabetes. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction Definitions and Classification of blood pressure levels according to 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension;

Category	Systolic (mmHg)		Diastolic (mmHg)
Optimal	<120	and	<80
Normal*3	120-129	and/or	80-84
High normal	130-139	and/or	85-89
Grade 1 hypertension	140-159	and/or	90-99
Grade 2 hypertension	160-179	and/or	100-109
Grade 3 hypertension	>180	and/or	>110
Isolated systolic hypertension	>140	and	<90

Isolated systolic hypertension should be graded (1, 2, 3) according to systolic blood pressure values in the ranges indicated, provided that diastolic values are <90mmHg

Source: The European Society of Hypertension and European Society of Cardiology Task Force Members. 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension. J Hypertens 2007; 25: 1751-1762

\*1 American Diabetes Association. The Diabetes-Heart Disease Link Surveying Altitudes Knowledge and Risk(2002)

## SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair by yourself and never try to disassemble the monitor under any circumstances.

### 1. Error messages

Message	Cause	What to do
Er1	No pulse signal detected	Check the wearing position of the cuff and remeasure it according to the required posture
Er2	Excessive interference	<ul style="list-style-type: none"> <li>Don't talk or shake your arms when measuring.</li> </ul> If error still remains, contact customer service for help.

Er3	Abnormal blood required posture pressure results	<ul style="list-style-type: none"> <li>• Check the wearing position of the cuff and remeasure it according to the required posture</li> </ul> If error still remains, contact customer service for help
ErP	Pressure is not up or too fast blocked	1. Check whether the air circuit is open or blocked 2. Check whether the cuff is too loose or too tight If error still remains, contact customer service for help
HI (Display in diastolic area)	Pressure over limit	Measure according to the required posture, do not squeeze the cuff during the measurement. If error still remains, contact customer service for help

## 2. Troubleshooting

1. If no display appears after pressing START button:

Possible cause	What to do
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed

2. If the heart rate is higher/lower than user's average:

Possible cause	What to do
Movement during measurement	Repeat measurement
Measurement taken just after exercise	Rest at least 30 minutes beforehand

3. If the result is higher/lower than user's average measurement:

Possible cause	What to do
May not be in correct position while measuring.	Adjust to the correct position to measure.
Blood pressure naturally varies from time to time	Keep in mind for next measurement.

4. If the cuff inflates again while measuring:

Possible cause	What to do
Cuff is not fastened.	Fasten the cuff again.
If user's blood pressure is higher than the pressure the device has inflated, the device will automatically increase the pressure and start to inflate again. Stay relaxed and wait for the measurement.	

## SYMBOL INFORMATION

Symbol	Referent
	Consult instructions for use
	Manufacturer
<b>SN</b>	Serial number
	Caution, consult accompanying documents
	Do not use if package is damaged
	Type B Equipment
<b>EC REP</b>	Authorized representative the European Community
	Dispose of in accordance with the requirements of your country.

## SYMBOL INFORMATION

Name	Full-automatic arm blood pressure monitor
Display mode	LCD Digital Display
Measurement method	Oscillographic measurement method
Measurement range	Pressure: 0~280mmHg (0kPa-37,3kPa); Pulse: 40~199 beats/min

Accuracy	Pressure: within +/-3 mmHg; Pulse: within +/-5%
Memory	2x99 (for two users each user 99 times memories)
Unit	mmHg/kPa
Power supply	USB DC 5V or 4 pcs AAA batteries; Input: AC 100-240V 50/60Hz; Output: DC 5V 1000mA
Automatic power off function	Without operation for 1 minute
Operating environment	5~40°C; 15%~80%RH
Storage temperature	-20~55°C; 10%~93%RH
Measuring range	About 220~320mm
Weight	About 430g
Dimensions	About 131x97x59mm
Security classification	Internal power supply type B
Pressure mode	Automatic pressure mode
Exhaust way	Automatic rapid exhaust
Pressure testing	Semiconductor pressure sensor
Pulse testing	Semiconductor pressure sensor

**Note:** If the above stored or used outside the specified temperature and humidity ranges, the system may not be able to achieve the above performance specifications.  
This device is not available for infant and it applies only to adults.  
This product has got the product registration certificate, if you have doubt about the effectiveness of product, you can search the company website to view certificate.

### Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

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